

easyBlood STARlet

The easyBlood STARlet assay ready workstation fully automates the process of blood fraction level detection and fraction transfer, creating the ideal starting sample for a wide range of assay applications that require whole blood components. On deck imaging of multiple blood tube types, coupled with pipetting, maximizes the number of tubes that can be imaged and processed at one time.



Popular Tube Types Supported on this Platform



BD™ Vacutainer® (13 x 75 mm, 4 mL;
13 x 100 mm, 6 mL; 16 x 100 mm, 10 mL)



Greiner Bio-One
VACUETTE® (6 mL, 9 mL)



Sarstedt® Monovette®
(4.9 mL, 7.5 mL)

Workflow

Blood Fraction Detection and Aliquoting

Autoload

Image

Analyze

Aliquot



hamiltoncompany.com/heretohelp

©2020 Hamilton Company. All rights reserved.
BD and Vacutainer are trademarks or registered trademarks of Becton, Dickinson and Company.
VACUETTE is a registered trademark of Greiner Bio-One GmbH.
Monovette® and Sarstedt® are registered trademarks of Sarstedt AG & Co. Nümbrecht, DE.
All other trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.
Lit. No. FL-2004-30 v1.0 – 05/2020